




SiDermique

Body Shape

Hyaluronic acid gel





SiDermique Body Shape is highly purified cross linked hyaluronic acid gel of excellent viscosity produced by a monophasic particle technology of triple level purification.

Its highest quality and excellent features are the best solution for modeling and correcting body contours and shapes.

→ Produced by a monophasic technology of triple purification

→ The concentration of purest hyaluronic acid is 24 mg/ml, which is the optimal concentration for injection into the middle and deep layers of the dermis, to form stable and beautiful shapes.

→ SiDermique Body Shape Gel is of non-animal origin, its high viscosity provides beautiful shapes and long-lasting effects with minimum swelling.

→ Safe and hypoallergenic.

→ Smooth and easy to use.

Indications for use



- restoration of volume and giving a beautiful shape to the parts of the body
- breast shape correction
- correction of skin relief after liposuction
- intimate contouring
- elimination of visible defects after surgeries or injuries

SIDE EFFECTS

Swelling and minor pain are the most common side effects experienced at the site of injection.

Rare perioperative infection within one breast is also another side-effect, with symptoms including tender pain and redness.

Antibiotics are usually prescribed to avoid this, and a pre-assessment should be carried out before the treatment to see if the patient has any existing skin infections.

SPECIFICATION



SiDermique Body Shape

Product type: Hard

Hyaluronic Acid gel of medical grade (Cross linked)

HA volume: 50ml / syringe

HA concentration: 20ml / ml

BDDE < 0.2ppm

Gel type: Mono-phasic

Package: 1 syringe / box

Store at a cool temperature at a cool dark place



SiDermique

THANK YOU

"You never know until you experience"